


ABSORBABLE SURGICAL SUTURE

TYPE	Twisted multifilament with a monofilament appearance	
COMPOSITION	Strands of purified collagen taken from the serosal layer of selected bovines, tanned with chrome salts. Origin from BSE-free classified countries.	
COATING	None	
TISSUE REACTION	Moderate	
ABSORPTION	By phagocytosis. The rate of absorption is slowed down by chromizing. Complete mass absorption in approximately 90 days. Approximately 50% of tensile strength remain after 14 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.	
TENSILE STRENGTH	Post implantation	Approximate % original strength remaining
	14 days	50 %
	28 days	20 %
SUTURE COLOUR	Dark brown	
RANGE	USP 6/0 - 4	
INDICATIONS	Catgut chrom sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, gastrointestinal surgery, gynaecology, obstetrics, urology, ophthalmic surgery	








Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	7	3	150 cm	brown		CC03WN150
	8	4	150 cm	brown		CC04WN150












Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR20 20 mm 	1.5	5/0	75 cm	brown		CC50HR20
	2	4/0	75 cm	brown		CC40HR20
	3	3/0	75 cm	brown		CC30HR20
HR26 26 mm 	2	4/0	75 cm	brown		CC40HR26
	3	3/0	75 cm	brown		CC30HR26
	3.5	2/0	75 cm	brown		CC20HR26
HR30 30 mm 	3.5	2/0	75 cm	brown		CC20HR30
HR40 40 mm 	3.5	2/0	75 cm	brown		CC20HR40
	3.5	2/0	90 cm	brown		CC20HR40-90
	4	0	75 cm	brown		CC00HR40
	4	0	90 cm	brown		CC00HR40-90
	5	1	75 cm	brown		CC01HR40
	5	1	90 cm	brown		CC01HR40-90
	6	2	75 cm	brown		CC02HR40
HR48 48 mm 	6	2	75 cm	brown		CC02HR48



Round bodied needle – taper point – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR30 30 mm 	4	0	75 cm	brown		CC00DR30-75
	4	0	90 cm	brown		CC00DR30
DR45 45 mm 	5	1	90 cm	brown		CC01DR45

 Reverse cutting needle – 3/8 circle							
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION
DS12	12 mm		1	6/0	75 cm	brown	CC60DS12-75
			1.5	5/0	75 cm	brown	CC50DS12-75
DS16	16 mm		2	4/0	75 cm	brown	CC40DS16
DS19	19 mm		1.5	5/0	75 cm	brown	CC50DS19
			2	4/0	75 cm	brown	CC40DS19
			3	3/0	75 cm	brown	CC30DS19
DS26	26 mm		3	3/0	75 cm	brown	CC30DS26
			3.5	2/0	75 cm	brown	CC20DS26
DS30	30 mm		3.5	2/0	75 cm	brown	CC20DS30
DS37	37 mm		4	0	75 cm	brown	CC00DS37