

LUXPET

NONABSORBABLE SURGICAL SUTURE
POLYESTER BRAIDED

TYPE	Braided multifilament
COMPOSITION	Polyester - a polymer of polyethylene terephthalate
COATING	Silicone
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue
SUTURE COLOUR	Green White
RANGE	USP 5/0 - 2
INDICATIONS	Luxpet sutures are intended for use in general soft tissue closing and/or ligation; especially in ophthalmic surgery, general surgery, skin closure, orthopaedics, gastrointestinal surgery, obstetrics



Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3.5	0	150 cm	green		PE00WN150
	4	1	150 cm	green		PE01WN150
	5	2	150 cm	green		PE02WN150

Without needles – ligapack

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	250 cm	green		PE30WN250
	3	2/0	250 cm	green		PE20WN250



Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16 16 mm 	1.5	4/0	75 cm	green		PE40HR16
HR20 20 mm 	2	3/0	75 cm	green		PE30HR20
HR26 26 mm 	2	3/0	75 cm	green		PE30HR26
	3	2/0	75 cm	green		PE20HR26
HR30 30 mm 	3	2/0	75 cm	green		PE20HR30
	3.5	0	75 cm	green		PE00HR30
HR36 36 mm 	3.5	0	75 cm	green		PE00HR36
	4	1	75 cm	green		PE01HR36
HR40 40 mm 	4	1	75 cm	green		PE01HR40
	5	2	75 cm	green		PE02HR40



Round bodied needle – taper point – 1/2 circle – double armed needle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR25 25 mm 	3	2/0	75 cm	green		PE20HR25-D

 Round bodied needle – tapercut point – 1/2 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm 	1.5	4/0	75 cm	green		PE40HRT18
	2	3/0	75 cm	green		PE30HRT18

 Round bodied needle – tapercut point – 1/2 circle – double armed needle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT17 17 mm 	3	2/0	90 cm	green		PE20HRT17-D

 Round bodied needle – tapercut point – 1/2 circle – double armed needle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT17 17 mm 	3	2/0	90 cm	white green	8 sutures / foil 5 foils / box with pledget 6 x 3 mm oval	PE20HRT17-DP
HRT25 25 mm 	3	2/0	90 cm	white green	8 sutures / foil 5 foils / box with pledget 6 x 3 mm oval	PE20HRT25-DP

 Reverse cutting needle – 3/8 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS16 16 mm 	1	5/0	75 cm	green		PE50DS16
	1.5	4/0	75 cm	green		PE40DS16
DS19 19 mm 	1.5	4/0	75 cm	green		PE40DS19
	2	3/0	75 cm	green		PE30DS19
DS24 24 mm 	2	3/0	75 cm	green		PE30DS24
	3	2/0	75 cm	green		PE20DS24
DS30 30 mm 	3.5	0	75 cm	green		PE00DS30