

STEEL

NONABSORBABLE SURGICAL SUTURE

STAINLESS STEEL MONOFILAMENT

TYPE	Monofilament
COMPOSITION	Iron, nickel and chromium alloy
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	Metal
RANGE	USP 1 - 7
INDICATIONS	Steel sutures are intended for use in closure of sternum, orthopaedics, abdominal wall closure and hernia repair



 Round bodied needle – tapercut point – 1/2 circle							
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT40	40 mm		6	3 & 4	4 x 45 cm	metal	ST03HRT40
HRT48	48 mm		7	5	4 x 45 cm	metal	ST05HRT48
			9	7	4 x 45 cm	metal	ST07HRT48

 Reverse cutting needle – 1/2 circle							
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS40	40 mm		4	1	4 x 45 cm	metal	ST01HS40
			5	2	4 x 45 cm	metal	ST02HS40
HS48	48 mm		6	3 & 4	4 x 45 cm	metal	ST03HS48
HS50	50 mm		7	5	4 x 45 cm	metal	ST05HS50
			8	6	4 x 45 cm	metal	ST06HS50
HS58	58 mm		7	5	4 x 45 cm	metal	ST05HS58
			8	6	4 x 45 cm	metal	ST06HS58

 Reverse cutting needle – 1/2 circle – PATELLA SET							
NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
NH120	120 mm		9	7	60 cm	metal	ST07HS120
						PATELLA SET	